

ALLURE

Total Oxygen Therapy Unit

You can have beauty inside and out with Allure - the world's first complete oxygen package offering a highly advanced beauticians' wand for oxygen facials and aromatherapy oxygen for inhalation!

The multi-purpose ALLURE – Total Oxygen Therapy Unit can be used in numerous locations for many purposes:

ALLURE OXYGEN:

Beauty Spas – for the ultimate in oxygen facials and hydrating tired skin

Health Clubs – great for energy and even better for de-stressing and relaxation

Airport Lounges – combat the affects of jet lag with revitalising oxygen

Home or Office – increase your brain activity, improve sleep patterns and general health

Many people do not know their bodies are suffering from low oxygen and toxins caused by the affects of pollution, bad diets, air travel, smoking, alcohol, infections, pollution, lack of exercise and low oxygen content in the air - but they may recognise the symptoms which include: low energy, interrupted sleeping patterns, frequent colds and infections, poor concentration and general poor health.

**RENTAL OR LEASING
OPTIONS AVAILABLE**

**IN AN EVER CHANGING WORLD
HEALTH IS EVERYTHING!**



The Allure - Total Oxygen Therapy Unit offers a new and unique way of combining oxygen inhalation with oxygen beauty facials. Delight your customers with two revolutionary treatments in one.

HOW IT WORKS:

The Allure - Total Oxygen Therapy Unit draws in ambient air containing 21% oxygen and increases the oxygen purity to around 90%. The unit is easy to operate and sits on wheels for easy portability!

BREATHING AROMATHERAPY OXYGEN:

With every breath your clients will inhale pure oxygen combined with sensual and revitalizing aromatherapy fragrances. Simply switch on, sit back and breathe for a totally relaxing and rejuvenating experience.

The benefits of oxygen are widely documented and include:

- Improved energy levels and vitality
- Clarity of the mind and increased concentration
- De-stressing and calming effects
- Relief from headaches and jet lag
- Detoxifying and cleansing
- Hydrating skin and reducing aging lines



OXYGEN FACIALS:

Delight your clients with the latest in facial therapy by simply connecting the beauty wand to the unit. The 'Allure' oxygen facial is the ultimate experience in non-surgical, non-invasive skin conditioning... combat those lines and let your skin breathe!

Oxygen treatment has received many accolades from physicians, beauty therapists, spa owners and celebrities alike. Oxygen facials have become so popular, even pop diva Madonna has been raving about the treatments in an interview with Harper's Bazaar. Kate Moss, Sienna Miller, Cameron Diaz and Claudia Schiffer are also reported devotees.

"It's the ultimate hydration therapy that makes your skin look better instantly and stay that way for days," said Dr. Bradford Katchen, a dermatologist in New York City, *"the facials [also] provide a moisture boost that makes skin smoother so that it is easier to apply make up."*

GENERAL SPECIFICATIONS:

- Quick Aroma Change Technology
- 8 Aroma Cartridge, 4 per person
- Supplied with Beauty Wand and Instructions
- Low operating costs
- Power capabilities for all power configurations (220v, 110v, 100v)
- No maintenance and no medical license required
- Telephone technical and set up helpline
- Marketing material supplied
- Nose hoses and connection kit
- Latest technology, whisper quiet oxygen concentrator
- One-year's manufacturers warranty
- CSA and CE certified



The aroma cartridges are designed to last a long time and come in many exciting fragrances. Here are a few of our most popular aromas:

Uplifting Wintergreen – stimulating and warming

Chillin Eucalyptus – purifying and energising

Invigorating Berry – cheering and cleansing

Purity Lemon Grass – revitalising and refreshing

The Beach wisteria – visualizing and exotic

Relax with juniper – nurturing and restoring

Zen Mint and Rosemary – cooling and sensual

Serenity Lavender – relaxing and calming

OXYGEN FACIAL TREATMENTS:

The Allure Total Oxygen Therapy Unit is user-friendly and safe. It is compatible with all skin types and uses oxygen that is 100% natural. An Allure treatment temporarily diminishes all the tiny imperfections that would otherwise be visible without any side effects. It is designed to be used alone or with skin care products of your choice. The pressurised oxygen speeds the skin's absorption of moisturising agents and produces impressive results. Achieve a dramatic reduction in the appearance of irritated, dry skin or just receive better circulation and a glow to tired, sallown skin in a specific area on the face, chest, back or any other area of the body.

ACNE TREATMENT:

Topical oxygen affects the epidermis and works to eliminate staphylococcus bacteria, which are partially responsible for causing acne pustules. The Allure Oxygen system wand pinpoints oxygen directly to specifically affected areas eliminating erythema (redness) and irritation. Oxygen is noted and documented as an anti-bacterial agent for infections and bacteria and also helps eliminate toxins trapped within the skin.

We provide oxygen solutions to an extensive global customer base in Europe, the Middle East, USA and Asia. Used for a variety of medical and recreational uses, our oxygen products benefit: hospitals; care homes; home users; health clubs; beauty salons; retailers; health organisations.

Our innovative products have improved the circumstances of all types of customers. Some of our customers have reported using oxygen to help with: **Sleep Apnea**; debilitating affects of **ME; MS; Allergies; Asthma; Sports**; affects of **Flying and Jet Lag**; Concentrating at **Work; General Health**.



The Allure beauty crème range



Pure O2 Ltd, The Old Chapel, Manchester Road, Carrington Village, Manchester, M31 4BL, UK
Website: www.healthoxygen.com Email: info@healthoxygen.com
UK Tel: 0870 712 02 02 INT Tel: 0044 161 776 2929

Disclaimer: Any practices or advice given or offered by Pure O2 is not intended to replace the services of your doctor, or to provide an alternative to professional medical treatment where it suggests the possible usefulness of certain practices in relation to certain illnesses or symptoms. It does so solely for educational purposes – either to explore the relationship of natural breathing to health, or to expose the viewer to alternative healing approaches. None of our statements are intended to diagnose, treat, cure or prevent disease. WARNING: Persons with any Medical Conditions or lung disorder including Asthma or Emphysema should consult their doctor before using recreational oxygen!