



Pulse On Demand

Continuous Flow

Battery Guide H=hours / M=minutes	FAA Approved	Weight kg	Trolley optional	Suitable for CPAP	1lpm	2lpm	3lpm	4lpm	5lpm	6lpm	0.5lpm	1lpm	1.5lpm	2lpm	3lpm
Airsep Focus with micro battery	YES	1.2	NO	NO	-	1h30m	-	-	-	-	-	-	-	-	-
Airsep Focus with extended battery	YES	1.9	NO	NO	-	3h30m	-	-	-	-	-	-	-	-	-
Airsep freestyle internal battery	YES	2	NO	NO	3h30m	2h30m	2h	-	-	-	-	-	-	-	-
Airsep Freestyle external battery	YES	2.7	NO	NO	10h	6h	3h	-	-	-	-	-	-	-	-
NEW Airsep freestyle clip on battery	YES	2	NO	NO	6h30m	3h30m	3h	-	-	-	-	-	-	-	-
NEW Airsep5 freestyle clip on battery	YES	2.7	NO	NO	4h15m	2h30m	2h	1h45m	1h15m	-	-	-	-	-	-
Inogen One G4 4 cell battery	YES	1.27	NO	NO	2h40m	2h15m	1h20m	-	-	-	-	-	-	-	-
Inogen One G4 8 cell battery	YES	1.5	NO	NO	5h	4h30m	2h30m	-	-	-	-	-	-	-	-
Inogen One G3 8 cell battery	YES	2.2	NO	NO	4h30m	4h	3h	2h	1h45m	-	-	-	-	-	-
Inogen One G3 16 cell battery	YES	2.7	NO	NO	9h30m	8h30m	6h	4h	3h20m	-	-	-	-	-	-
Philips Simply Go MINI 8 cell battery	YES	2.3	NO	NO	4h30m	4h30m	2h30m	2h	1h30m	-	-	-	-	-	-
Philips Simply Go MINI 16 cell battery	YES	2.7	NO	NO	9h	9h	5h	4h	3h	-	-	-	-	-	-
Philips Simply Go battery	YES	4.5	YES	YES	3h40m	3h30m	3h	2h40m	2h15m	1h55m	3h	2h15m	1h30m	0h50m	-
DeVilbiss i-Go battery	YES	10	YES	YES	5h25m	4h40m	4h	3h30m	3h10m	3h	-	4h	-	2h25m	1h35m
Gce Zeno	YES	4.66	YES	YES	4h	4h	3h	2h15m	2h	1h45m	3h	1h45m	1h15m	0h45m	-
Sequal Eclipse battery	YES	8.35	YES	YES	5h25m	5h	4h50m	4h	3h40m	3h30m	4h	3h35m	2h30m	2h	1h25m
Sequal Equinox with 12 cell battery	YES	6.35	YES	YES	2h30m	2h30m	2h15m	2h15m	1h30m	1h30m	2h15m	2h15m	1h45m	1h15m	0h45m
Sequal Equinox with 24 cell battery	YES	7.72	YES	YES	5h45m	5h45m	5h15m	5h15m	3h45m	3h45m	4h45m	4h45m	3h30m	2h45m	1h45m